



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



STATIM 5000S autoclave (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Made in Canada,
Power supply: 230 V,
Frequency: 50/60 Hz,
STATIM 5000s
is a handheld autoclave for hospitals and medical offices of many specialties, ideal for the rapid sterilization of ophthalmic, surgical, endodontic, and other delicate instruments,
Class: S,
Operating noise: 56 dB,
Features:
Energy-efficient,
Fast sterilization cycles - up to 9 minutes for unwrapped instruments,
Low distilled water consumption,
Simple design for easy service and maintenance,
Compliance with the EN13060 standard for sterilization of type A and B loads,
Sterilization programs:
Cycle for wrapped 134°C - 18°C min,
Unwrapped cycle 134°C - 9 min,
Plastics and rubber 121°C - 24 min,
Drying - minimum 60 min,
Sterilization and drying of wrapped products - approximately 82 min,
Sterilization and drying of unwrapped products - approximately 45 min,
Menu languages: PL, ENG, GER,
Includes:
STATIM cassette,
Condensate container,
New biological filter,
Instruction manual in English (PDF),
Power cord,
Dimensions:
Autoclave: 41 x 60 x 19 cm,
Internal cassettes: 18 x 37 x 6.5 cm,
External cassettes: 20 x 50 x 8 cm,
Weight: 34 kg,
Cat 02,
Valid inspection and ready for operation,
Technical Passport issued (Service Report) valid for 12 months,
Warranty:
6 months for the domestic market (Poland),
3 months for the international market,
Warranty can be extended for an additional fee to 12 months for the domestic and international markets,
Financing options (Poland only): Installments, Leasing, Loan,

