



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **STATIM 5000S ENDO autoclave (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Actual product photos,  
Made in Canada,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
STATIM 5000s  
is a handheld autoclave for hospitals and medical offices of many specialties, ideal for the rapid sterilization of ophthalmic, surgical, endodontic, and other delicate instruments,  
Class: S,  
Operating noise: 56 dB,  
Features:  
Energy-efficient,  
Fast sterilization cycles - up to 9 minutes for unwrapped instruments,  
Low distilled water consumption,  
Simple design for easy service and maintenance,  
Compliance with the EN13060 standard for sterilization of type A and B loads,  
Sterilization programs:  
Cycle for wrapped 134°C - 18°C min,  
Cycle for unwrapped items 134°C - 9 min,  
Plastics and rubber 121°C - 24 min,  
Drying - minimum 60 min,  
Sterilization and drying of packaged items - approximately 82 min,  
Sterilization and drying of unwrapped items - approximately 45 min,  
Menu in: Polish,  
Includes:  
ENDO STATIM cassette,  
Condensate container,  
New biological filter,  
Instruction manual in English (PDF),  
Power cord,  
Dimensions:  
Autoclave: 41 x 60 x 19 cm,  
Internal cassettes: 18 x 38 x 7.5 cm,  
External cassettes: 19.5 x 49.5 x 8 cm,  
Weight: 34 kg,  
Cat 03,  
Valid inspection and ready for operation,  
Technical Passport issued (Service Report) valid for 12 months,  
Warranty:  
6 months for the domestic market (Poland),  
3 months for the international market,  
Warranty can be extended for an additional fee to 12 months for the domestic and international markets,  
Financing options (Poland only): Installments, Leasing, Loan

