



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Shock wave (Reconditioned)**

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Actual product photos,

Manufactured by P.R.C.,

Power supply: 230 V,

Frequency: 50/60 Hz,

The therapy uses acoustic sound waves generated inside the head and gently transferred to the body areas affected by cellulite. By using shockwaves to break down subcutaneous fat tissue, this treatment is an excellent alternative to traditional methods.

Applications:

Cellulite treatment,

Localized obesity removal,

Removal of loose skin,

Removal of skin irregularities after liposuction,

Reduction of limb swelling,

Reduction of stretch marks,

Technical specifications:

Technology: shockwave,

Power/energy level: 200 mj,

Frequency: 1 - 16 Hz,

2 operating modes,

Menu language selection: English, Spanish,

Includes:

Power cord,

Handle with head,

Dimensions: 31.5 x 28 x 18 cm,

Weight: 4.2 kg,

Updated inspection and ready to use,

Issued Technical Passport (Service Report) valid for 12-month warranty

Warranty:

6 months for the domestic market (Poland),

3 months for the international market

Possibility of extending the warranty for an additional fee to 12 months for the domestic and international markets

Financing options (Poland only): Installments, Leasing, Loan

