



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **OLYMPUS Evis Exera II CV-180 Video Processor with CLV-180 Light Source (Kit) (Reconditioned)**

Used,

Technical condition: very good,

Visual condition: very good,

Japanese production,

The set includes:

OLYMPUS EVIS EXERA II Model CV-180 image processor:

Power supply: 230 V,

Frequency: 50/60 Hz,

Power: 150 VA,

The unit has the ability to process high-definition HDTV images to provide the best image quality in endoscopy and laparoscopy,

Compatible with EVIS 100/130/140 series and EVIS EXERA 160 series endoscopes,

HD / SD SDI output for high quality video transfer,

Has a port for xD memory cards,

Foiled control panel,

User manual in English (PDF),

Dimensions: 50 x 37 x 9 cm,

Weight: 10 kg,

OLYMPUS EVIS EXERA II Light Source Model CLV-180:

Power supply: 230 V,

Frequency: 50/60 Hz,

Power: 500 VA,

Adapted for flexible endoscopes,

Main lighting: xenon bulb MD 631 / 300 W,

Emergency lighting: 12 V / 35 W halogen bulb,

Foiled control panel,

User manual in English (PDF),

The unit has minor signs of use in the form of a dent in the cover,

Dimensions: 54 x 39 x 16 cm,

Weight: 15 kg,

OLYMPUS MAJ-1428 Keyboard,

Light Control Cable MAJ-1411,

Power cable - 2 pcs,

Pigtail MAJ-1430 cable not included - can be retrofitted at additional cost,

Weight of the set: 27 kg,

Cat 03,

Has a current Technical Passport issued,

Warranty:

6 months for domestic market - Poland,

3 months for the international market,

Possibility of financing (Poland only): Installments, Leasing, Loan,

