



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **HAIER BIO MEDICAL DW-86.W100 low-temperature freezer (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Actual product photos,  
Manufactured by P.R.C.,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 650 W,  
Freezers designed for long-term storage of biological material, e.g., viruses, erythrocytes, leukocytes, tissues, and laboratory reagents,  
Wide tolerance for supply voltage values - operation possible even with unstable power supplies,  
Wide tolerance for ambient temperature; ambient temperature range required for proper device operation: +10°C to +32°C,  
Specifications:  
Capacity: 100 l,  
Operating temperature range: -40°C to -86°C,  
Temperature hysteresis (temperature deviation from actual temperature): 5°C,  
Controller type: Microprocessor-based,  
Display: LED,  
Temperature programming from -40°C to -86°C,  
Programming alarm thresholds for too low/too high temperatures,  
Temperature alarm for too low/too high temperatures,  
Power failure alarm,  
Temperature sensor damage alarm,  
Power supply voltage alarm for too high/too low temperatures,  
Ambient temperature alarm for too high temperatures,  
Condenser temperature alarm for too high temperatures,  
Low battery charge alarm for the control panel during power outages,  
Includes:  
Power cord,  
Tiered basket for storing test tubes,  
7 (pcs.) test tube storage packages,  
Dimensions:  
External: 77 x 82.5 x 111.2 cm,  
Internal: 46 x 45 x 50 cm,  
Weight: 138 kg,  
Currently inspected and ready for operation,  
Issued Technical Passport (Service Report) valid for a period of 12-month warranty  
Warranty:  
6 months for the domestic market (Poland),  
3 months for the international market  
Possibility of extending the warranty for an additional fee to 12 months for the domestic and international markets  
Financing options (Poland only): Installments, Leasing, Loan

