



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **OLYMPUS UCR Insufflator (Reconditioned)**

Used,  
Technical condition: very good,  
Visual condition: very good,  
Made in Japan, (2012)  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 40 VA,  
The OLYMPUS UCR insufflator  
increases patient comfort during and after colonoscopy,  
Carbon dioxide is absorbed by the body 150 times faster than air,  
Easy operation with a single button,  
Pressure level display,  
Three controlled airflow settings available,  
Max supply pressure: 45 kPa,  
Can be retrofitted for an additional fee with the following accessories shown in the last photo:  
Water bottle (Olympus MAJ-902),  
Gas tube (Olympus MAJ-1741),  
Includes:  
CO2 cylinder hose,  
Power cord,  
Dimensions: 12.5 x 30 x 15 cm,  
Weight: 5 kg,  
Valid Technical Passport issued,  
Warranty:  
6 months for the domestic market - Poland,  
3 months for the international market,  
Financing options (Poland only): Installments, Leasing, Loan,

