



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



MINDRAY IPM12 Patient Monitor/Cardiac Monitor (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Made in Germany,
Power supply: 230 V,
Battery backup: via two built-in 11.1 V / 4.5 Ah batteries (1 new one included) -
(possibility of retrofitting a second new one for an additional fee),
Frequency: 50/60 Hz,
Technical specifications:
12.1" color TFT LCD screen,
Compact and modular design - measurement functions can be expanded by adding modules,
Screen layout configuration: order, number, and color of curves; display of short trends next to dynamic curves; oxycardiorespirogram screen;
large digit screen; user-configurable and memorizable custom screens,
Extensive data storage and review functions: 120-hour graphical and tabular trends, 1000 non-invasive blood pressure measurement results,
100 alarm events, 100 arrhythmia events, up to 48 hours of real-time waveform storage (full disclosure),
Alarm system enabling audible and visual alarm notification; individual, manual alarm settings for individual parameters and automatic setting
of alarm limits based on current parameter values,
Printing function for data, waveforms, and trends on a standard laser printer on A4 paper,
Built-in thermal printer with a 5 cm paper width,
Measurement parameters:
ECG,
SPO2,
NIBP,
Temperature [T1 and T2],
Patient monitor equipped with a CO2 gas measurement module - no water trap included (can be retrofitted for an additional fee),
Menu language: Polish,
Includes:
User's manual in English (PDF),
Power cord,
New ECG cable 5-lead,
New NIBP air cuff with airway,
New SPO2 clip sensor,
New connecting cable,
Dimensions: 30 x 15 x 30 cm,
Weight: 8 kg,
Updated inspection and ready for use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6 months for the domestic market (Poland),
3 months for the international market,
Possibility of extending the warranty for an additional fee to 12 months for the domestic and international markets,
Financing options (Poland only): Installments, Leasing, Loan,

