



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



CARBOLITE PIF60 laboratory incubator (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Made in Germany,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 675 W,
The CARBOLITE PIF60
is a laboratory device with natural air circulation. They provide basic drying conditions, such as long-lasting, precise high temperature control.
Specifications:
Natural air circulation,
Digital temperature setting up to 80°C (accurate to 1°C),
Air circulation control using a knob,
Volume: 66 liters,
Includes:
Power cord,
(1 pc.) tray,
Dimensions:
Internal: 39 x 42 x 40 cm,
External: 76.5 x 57 x 57 cm,
Weight: 45 kg,
Updated inspection and ready to use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6 months for the domestic market (Poland),
3 months for the international market,
Possibility to extend the warranty for an additional fee to 12 months for the international market. domestic and international,
Financing options (Poland only): Installments, Leasing, Loan,

