



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Treatment lamp MAQUET LUCEA 10 (ceiling) (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Made in Sweden,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 22 VA,
Maquet Lucea 10
is a ceiling-mounted treatment lamp.

Features:
Thanks to its excellent color reproduction, it is particularly suitable for dermatology.
The use of LEDs prevents the light from emitting heat, which helps distinguish healthy veins from varicose veins.
The lampshade can be adjusted to any position thanks to multi-level adjustment.
It has a two-part arm with a total length of 80 cm.
The lampshade is equipped with 10 LEDs that generate a maximum light intensity of 50,000 lx.

Specifications:
Color temperature: approx. 4500 K.
Color Rendering Index: 95 Ra,
Light field: Ø 13 cm,
Includes:
Lamp suspension column,
Masking frame,
Mounting kit,
Power supply,
Dimensions:
Shade: 22 x 17 cm,
Overall dimensions: 72 x 26 x 15 cm,
Length of the lamp suspension column: 100 cm,
Total weight: 11.5 kg,
Updated inspection and ready for operation,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6 months for the domestic market (Poland),
3 months for the international market,
Possibility of extending the warranty for an additional fee to 12 months for the domestic and international markets,
The warranty period does not cover bulbs and other light sources installed in the device,
Financing options (Poland only): Installments, Leasing, Loan,

