



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Polometer (perimeter) CARL ZEISS HUMPHREY 750 (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Made in the U.S.A.,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 400 W,
Halogen bulb: 12 V / 20 W,
Specifications:
Examination type: static and kinetic,
Field of view: 90°,
Maximum stimulus intensity: 10,000 ASB,
Cover illumination: 31.5 ASB,
Pulse duration: 200 ms,
Examination distance: 30 cm,
Equipped with a built-in 12" touchscreen monitor,
Electrical chin rest adjustment: up, down and sideways,
Hard drive and floppy disk drives,
Menu language selection: ENG, GER, ESP, FR, IT, JP, PT, SWE,
Included:
User's manual in English (PDF),
Power cord (2 pcs),
Printer cable,
LEXMARK MS312DN laser printer,
Patient button,
Dimensions:
Perimeter: 47 x 52 x 60 cm,
Printer: 40 x 38 x 28 cm,
Weight:
Perimeter: 36.4 kg,
Printer: 14 kg,
Updated inspection and ready to use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6 months for the domestic market (Poland),
3 months for the international market,
Possibility to extend the warranty for an additional fee to 12 months for the domestic market and international,
The warranty period does not cover light bulbs and other light sources installed in the device,
Financing options (Poland only): Installments, Leasing, Loan,

