



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Oxygen (O2) reducer VTI SELECT ERGO STANDARD 25 (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- German production,
- **Specification:**
  - Housing and fasteners: CuZn alloy / brass
  - Piston tightness: PA 6
  - O-rings: silicone, EPDM
  - Filter: Sintered Bronze
  - 12-step adjustable setting 0-25 l/min: 0 / 0.5 / 1 / 1.5 / 2 / 3 / 4 / 6 / 8 / 12 / 15 / 25,
  - 9/16" outlet with threaded hose connector,
  - Pressure: up to 200 bar,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
  - **6 months for the domestic market (Poland),**
  - **3 months to the international market,**
  - **Possibility of extending the warranty for an additional fee up to 12 months for the domestic and international markets,**
- Financing options (Poland only): Installments, Leasing, Loan,

