



**PRAISTON Sp. z o.o.**  
ul. Górska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67  
E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



## **EMED ES350 ARGON surgical diathermy (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Actual product photos,  
Made in Poland,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 500 W,  
Features:  
Generator operating frequency: 333 kHz,  
Surgical diathermy integrated with an Argon module,  
Allows operation in standard monopolar and bipolar cutting and coagulation modes, as well as in argon cutting and coagulation modes,  
Excellent in both open surgery and endoscopic and laparoscopic procedures,  
Laminated, touchscreen, and easy-to-read control panel,  
It has the following functions:  
**MONOPOLAR CUTTING:**  
Pure cutting 400 [W] for 200 [Ω],  
Mixed cutting I 180 [W] for 200 [Ω],  
Mixed Cutting II 150 [W] for 200 [Ω],  
Mixed Cutting III 150 [W] for 200 [Ω],  
Urological Cutting 9 Levels,  
Argon Shielded Cutting 350 [W] for 200 [Ω],  
Endoscopic Cutting 9 Levels,  
**MONOPOLAR COAGULATION:**  
Soft Coagulation 180 [W] for 50 [Ω],  
Forced Coagulation 180 [W] for 300 [Ω],  
Spray Coagulation 80 [W] for 1250 [Ω],  
Hybrid Coagulation 180 [W] for 200 [Ω],  
Argon Coagulation 80 [W] for 1250 [Ω],  
Pulsed Argon Coagulation 80 [W] for 1250 [Ω],  
**BIPOLAR CUTTING:**  
Bipolar cutting 150 [W] for 200 [Ω]  
(4 hemostasis levels),  
Urological cutting 9 levels,  
**BIPOLAR COAGULATION:**  
Bipolar coagulation 120 [W] for 50 [Ω],  
AutoStart / AutoStop automatic start and stop,  
**ARGON:**  
Argon type: pure argon 4.8 (99.998 [%]) or higher,  
Input pressure: 0.3-0.5 [MPa],  
Gas flow: 0.1-9.9 [l/min], adjustable in 0.1 [l/min] increments throughout the entire range,  
Safety systems:  
PowerStart system - supports initiation of the cutting process,  
SURGILOGIC system - automatic power adjustment to the tissue,  
NEM system - allows for contact quality control neutral electrode with the patient's skin, which helps eliminate burns,  
AUTOTEST - device and accessory test after powering on, displaying service codes,  
Includes:  
Bipolar cable for BOWA forceps  
(can be retrofitted with forceps and bipolar paste for an additional fee),  
Brand new monopolar electrode cable with BOWA holder - two buttons, 4 mm shaft  
(electrodes not included, can be retrofitted for an additional fee),  
Brand new BOWA clip cable for disposable neutral electrodes  
(electrodes not included, can be retrofitted for an additional fee),  
New cable foot switch,  
Power cord,  
Can be retrofitted with cables and argon probes for an additional fee,  
Dimensions: 41.5 x 42 x 18 cm,  
Weight: 10 kg,  
Currently inspected and ready to use,  
Technical Passport issued (Report Service) valid for 12 months,  
Warranty:  
6 months for the domestic market (Poland),  
3 months for the international market,  
Warranty can be extended for an additional fee to 12 months for the domestic and international markets,  
Financing options (Poland only): Installments, Leasing, Loan,

