



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MEDITOP MT750 ENT unit (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- Dutch production,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 2400 VA,
- **The MT750** is a mobile and compact ENT unit. The MT750 offers various storage options for instruments and supplies for minor procedures. It's an ideal treatment module for use in multi-functional spaces.
- The MT750 is equipped with 3 parking handles as standard, including for suction, compressed air and cold light,
- Freely extending drawers for accessories,
- Built-in surgical suction,
- **Unit equipped with:**
 - **NEOVO X-19AV Monitor:**
 - Screen size: 19" inches,
 - Screen format: 4:3,
 - Resolution: 1280 x 1024 px.,
 - Static/dynamic contrast: 700:1,
 - Screen brightness: 300 cd/m2,
 - Matrix response time: 4 ms,
 - Vertical/horizontal viewing angle: 150°/135°,
 - Connectors: DVI x 1 / VGA x 1,
 - **Surgical diathermy ERBE ICC 50 INT:**
 - Reconditioned (used),
 - Technical condition: very good,
 - Visual condition: very good,
 - Real product photos,
 - German production,
 - Power supply: 230 V,
 - Frequency: 50/60 Hz,
 - Power: 50 W,
 - **Technical specifications:**
 - Digital displays for cutting and coagulation power,
 - Foiled, easy-to-use control panel,
 - **It has the following functions:**
 - Cut (CUT) - 50 W at 500 Ω:
 - EFFECT 1 (low coagulation effect),
 - EFFECT 2 (high coagulation effect),
 - Coagulation (COAG) - 50 W at 200 Ω:
 - Forced (FORCED),
 - Soft (SOFT),
 - Bipolar (BIPOLAR),
 - **MEDITOP MT LUMITOP 150 light source,**
 - Reconditioned (used),
 - Technical condition: very good,
 - Visual condition: very good,
 - Real product photos,
 - Dutch production,
 - Power supply: 230 V,
 - Frequency: 50/60 Hz,
 - Power: 175 VA,
 - **OLYMPUS OTV-SX2C endoscopic camera image processor:**
 - Reconditioned (used),
 - Technical condition: very good,
 - Visual condition: very good,
 - Real product photos,
 - Japanese production,
 - Power supply: 230 V,
 - Frequency: 50/60 Hz,
 - Power: 50 VA,
 - Bulb: 15V / 150W (2 pcs.),
 - **Specification:**
 - TV type: PAL,
 - Color saturation adjustment: red, blue,
 - WHITE BALANCE function,
 - Aperture adjustment,
 - Video Enhancement,
 - **Outputs:**
 - Video OUT,
 - 2 x Y/C,
 - RGB,
 - H.COPY,
 - SYSTEM,
 - VTR RENOVATION,
 - **OLYMPUS OTV-SX2 camera head adapted to rigid endoscopes, placed in a portable case,**
- **Included in the set:**
 - ERBE ICC 50 INT user manual in English (PDF),
 - Power cord,
 - Fiber optic cable to the light source,
 - SERRES 2 l autoclavable jar with disposable insert,
 - OLYMPUS OTV-SX2 camera head with dedicated carrying case,
 - ERBE wired foot manipulator,
- Overall dimensions: 100 x 60 x 155 cm,
- Total weight: 120 kg,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months to the international market,**
 - **Possibility of extending the warranty for an additional fee up to 12 months for the domestic and international markets,**
 - **The warranty period does not cover light bulbs and other light sources installed in the device,**
- Financing options (Poland only): Installments, Leasing, Loan,

