



MELAG VACUKLAV 24 autoclave (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Made in Germany,
Power supply: 230 V,
Frequency: 50/60 Hz,
Power: 2200 W,
Specifications:
Temperature range up to 144°C,
Pressure up to 3 bar,
Operating time adjustable with a knob within a 70-minute range,
Autoclave capacity: 18.7 l,
2 operating modes:
20 minutes – 3 bar,
45 minutes – 1 bar,
Includes:
Power cord,
Tray basket,
2x 18 x 41 cm trays,
Tray holder,
Overall dimensions: 45 x 60 x 45 cm,
Weight: 30 kg,
Currently inspected and is ready for operation.
Issued Technical Passport (Service Report) valid for 12 months.
Warranty:
6 months for the domestic market (Poland),
3 months for the international market.
Possibility of extending the warranty for an additional fee to 12 months for the domestic and international markets.
Financing options (Poland only): Installments, Leasing, Loan.

