



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **WELCH ALLYN PROPAQ CS 242 patient monitor/cardiac monitor (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real photos of the product,  
Manufactured by U.S.A.,  
Power supply: 230 V,  
Emergency power supply: through 2 new built-in 8 V / 2.7 Ah batteries,  
Frequency: 50/60 Hz,  
PROPAQ CS 242  
is a compact patient monitor that allows you to control many of the patient's vital signs even under extreme conditions,  
The device has:  
8.5" LCD touchscreen display,  
The ability to connect to a defibrillator via the DEFIB SYNCHRO jack,  
Audible and visual alarms,  
Handle for moving and hanging,  
Equipped with the following modules with full wiring and sensors:  
ECG / ECG RESP - New 3 - lead cable with snap type connectors,  
SOP  
2,  
NIBP - New cuff with air line,  
TEMP 1/2 - No sensor included  
(can be retrofitted at additional cost),  
English language menu,  
Included:  
User manual in English (PDF),  
Power cord,  
19 V power supply,  
Dimensions: 24 x 20 x 20 cm,  
Weight: 6 kg,  
Cat 03,  
Has a current inspection and is ready for operation,  
Issued Technical Passport (Service Report) valid for 12 months,  
Warranty:  
6-months for domestic market (Poland),  
3-months for the international market,  
Possibility to extend the warranty up to 12-months for domestic and international market for a surcharge,  
Financing options (Poland only): Installments, Leasing, Loan,

