



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



WELCH ALLYN PROPAQ CS 244 Patient Monitor/Cardiac Monitor (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Made in the U.S.A.,
Power supply: 230 V,
Battery power supply: via two new built-in 8 V / 5 Ah batteries,
PROPAQ CS 242
is a compact patient monitor that allows you to monitor many of the patient's vital signs even in extreme conditions.
The device features:
An 8.5" LCD touchscreen display,
Connection to a defibrillator via the DEFIB SYNCHRO socket,
Acoustic and visual alarms,
Handle for carrying and hanging,
Equipped with the following modules with full wiring and sensors:
ECG / EKG RESP - 3-lead cable,
New SPO2 sensor,
NIBP - New cuff with air tube,
TEMP 1/2 - No sensors in the Included
(can be retrofitted for an additional fee),
Includes:
User's manual in English (PDF),
Monitor power supply,
3-lead ECG cable,
New SPO2 clip sensor,
New NIBP cuff with airway,
Dimensions: 24 x 20 x 20 cm,
Weight: 6 kg,
Updated inspection and ready for use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6 months for the domestic market (Poland),
3 months for the international market,
Possibility of extending the warranty for an additional fee to 12 months for the domestic and international markets,
Financing options (Poland only): Installments, Leasing, Loan,

