



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



CARL ZEISS KSK 150-FC colposcope (video track) (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- German production,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 150 VA,
- **ZEISS KSK 150-FC** is a mobile colposcope equipped with a vision track,
- Light source: 12V / 100W halogen bulb,
- Smooth regulation of light intensity,
- Fiber optic,
- Green filter,
- Magnification: 0.4x / 0.6x / 1x / 1.6x / 2.5x,
- Focal length (f): 300,
- Eyepiece: 10x,
- **Smooth regulation:**
 - Eyepiece heights in the range of 74 -120 cm,
 - Adjusting the distance between the glasses using the knobs,
- There are knobs on the pantograph arm that lock it in the selected position,
- **Colposcope extended with a vision track, which includes:**
 - **Microscope camera:**
 - Sensor: 48 million pixels CMOS,
 - Interface: HDMI (Type A). USB (type B).,
 - HDMI image resolution: 48 MP,
 - Video Resolution: 2K at 30fps 1080p at 60fps,
 - HDMI Resolution: 1920x1080p 60fps,
 - HDMI Scale: Support 4:3/16:9/16:10 screen,
 - USB 2.0 (driveless): 1920 x 1080 @ 30 fps 1280 x 720 @ 30 fps 640 x 480 @ 30 fps
 - **MSI PRO273U Monitor:**
 - Screen diagonal: 27" inches,
 - Matrix coating: matte,
 - Screen backlight: LED,
 - Matrix type: IPS,
 - Screen resolution: 3840 x 2160 px.,
 - Screen aspect ratio: 16:9,
 - Image refresh rate: 60 Hz,
 - Color space reproduction: Adobe RGB: 93% / DCI-P3: 98% / sRGB: 130%,
 - Number of displayed colors: 1.07 billion,
 - HDR: HDR Ready,
 - Matrix response time: 4 ms,
 - **CARL ZEISS light source:**
 - Light source: halogen lamp 12 V / 150 W,
 - Possibility of smooth regulation of light intensity,
- The structure operates on a mobile base equipped with 5 wheels (2 with brakes),
- **Included in the set:**
 - KSK FC-150 colposcope user manual in English (PDF),
 - Power cord,
- Dimensions: 135 x 110 x 70 cm,
- Weight: 68 KG,
- Cat 10,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months to the international market,**
 - **Possibility of extending the warranty for an additional fee up to 12 months for the domestic and international markets,**
 - **The warranty period does not cover light bulbs and other light sources installed in the device,**
- Financing options (Poland only): Installments, Leasing, Loan,

