



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **MEDTRONIC BIPHASIC LIFEPAK 12 Defibrillator (Cat. 81) (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: good,  
Actual product photos,  
Made in the U.S.A.,  
Power supply: via replaceable 12V 1.2 Ah batteries,  
(one included, second can be added for an additional fee)  
Frequency: 50/60 Hz,  
Power: 180 W,  
Menu language selection: EN, FR, DE, ES, SW, IT, FIN, NL, PL, DK, NR,  
Control panel keypad in: English,  
Built-in printer with paper width: 11 cm,  
Equipped with a carrying handle,  
Functions:  
SYNC - Synchronization - Cardioversion,  
ADVISORY - Activates the Continuous Patient Surveillance System (CPSS),  
ANALYZE - Activates the advisory system SAS (Shock Advisory System),  
Equipped with the following modules with wiring and sensors:  
ECG - 10-lead cable,  
SPO2,  
Cuff with tubing,  
AED electrodes,  
AED cable,  
For an additional fee, the defibrillator can be equipped with:  
Hard paddles,  
Additional battery,  
External charger from AD Elektronik,  
Includes:  
Instruction manual in Polish (PDF),  
Additional transport bags for accessories,  
Dimensions:  
Defibrillator: 20 x 40 x 34.5 cm,  
With bag: 23 x 62 x 34.5 cm,  
Weight: 8.5 kg,  
Cat 81,  
Certified and ready for use,  
Technical Passport (Report issued Service) valid for 12 months,  
Warranty:  
6 months for the domestic market (Poland),  
3 months for the international market,  
Warranty can be extended for an additional fee to 12 months for the domestic and international markets,  
Financing options (Poland only): Installments, Leasing, Loan,

