



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MEDTRONIC LIFEPAK 20 Defibrillator (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: good,
- **Real product photos,**
- Made in the USA,
- Power supply: 230 V,
- Emergency power supply: via 1 new built-in battery,
- Frequency: 50/60 Hz,
- Power: 120 VA,
- **Possibility to select menu language:** PL, EN, FR, DE, ES, SV, IT, DA, PT, ZH, CS, RU, JA, HU, KO, NL, NO, FI ,
- **Lifepak 20** is a biphasic, portable defibrillator that provides a complete response system for acute cardiac emergencies.
- **Technical specifications:**
 - **AED mode,**
 - **Synchronization,**
 - **Monitoring ECG parameters,**
 - Easy-to-read 6.2" display
 - Built-in printer (paper width: 50 mm),
- **Included in the set:**
 - User manual in English (PDF),
 - 3-lead ECG cable,
 - QUICK-COMBO cable,
 - Power cord,
- Dimensions: 26 x 24 x 21 cm,
- Weight: 6.5 kg,
- Cat D152,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months to the international market,**
 - **Possibility of extending the warranty for an additional fee up to 12 months for the domestic and international markets,**
- Financing options (Poland only): Installments, Leasing, Loan,

