



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



PHILIPS INTELLIVUE MP50 ANESTHESIA patient monitor/cardiomonitor (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Made in Germany (2012),
Power supply: 230 V,
Backup power supply: via 2 built-in batteries – new batteries can be retrofitted for an additional fee,
Frequency: 50/60 Hz,
The IntelliVue MP50 ANESTHESIA
is a patient monitor that monitors multiple vital signs.
Features:
A 12-inch TFT LCD panel display with SVGA resolution.
Minor scratches that do not affect the monitor's operation.
Operation via a touchscreen and a navigation joystick allows direct access to many functions using simple commands. The screen can display up to 6 waveforms, while the 12-lead ECG screen can display 12 ECG waveforms,
A comfortable handle and durable housing make it easy to carry,
Custom display options allow you to view and analyze data in graphical and numerical formats, collate measurements and trends in real time, and freely configure every screen element - from waveforms to parameter names,
IntelliVue technology - access information from various systems in the hospital network - pharmacy, laboratory, ECG management, radiology PACS, server-based hospital systems, and others - and deliver it to the patient monitor,
No separate hard drive or fan - IntelliVue monitors are small, quiet, stable, and start up quickly,
The navigation knob makes it easy to enter information and navigate the screen,
Audible and visual alarms,
Defib Sync port, LAN,
Equipped with:
PHILIPS M1029A Temperature Module:
The Philips temperature module provides continuous body temperature readings in various areas,
It allows you to monitor temperatures ranging from -1°C to 45°C with alarm thresholds that can be set to 0.1°C increments.
PHILIPS M1006B Invasive Blood Pressure Module:
Provides waveforms and numerical values and works with a wide range of pressure catheters and transducers.
Invasive pressures can be collected, via the module or an internal function, from adult, pediatric, and neonatal patients for display on IntelliVue patient monitors.
The Invasive Blood Pressure Module provides real-time waveforms and numerical values for systolic, diastolic, and mean blood pressure, as well as associated pulse rates.
PHILIPS IntelliVue X2 Measurement Module:
ECG/RESP – Brand new 3-lead cable.
SPO2 – Oxygen saturation sensor.
NBP – Cuff with air tube.
TEMP – Temperature connector port.
(sensor not included, can be retrofitted for an additional (subject to a fee),
Menu in: Polish,
Includes:
Instruction manual in Polish (PDF),
Power cord,
New 3-lead ECG cable,
SPO2 sensor,
NIBP cuff with air tube,
Dimensions: 22.5 x 36.5 x 34 cm,
Weight: 10 kg,
Updated inspection and ready to use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6 months for the domestic market (Poland),
3 months for the international market,
Possibility of extending the warranty for an additional fee to 12 months for the domestic and international markets,
Financing options (Poland only): Installments, Leasing, Loan,

