



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



WELCH ALLYN PROPAQ ENCORE 206 EL patient monitor/cardiac monitor (Reconditioned)

Used,
Technical condition: very good,
Visual condition: very good,
Manufacture U.S.A.,
Power supply: 230 V and through 2 new built-in 8 V / 2.7 Ah batteries,
Frequency: 50/60 Hz,
The Propaq Encore patient monitor
has been designed specifically for ambulance and critical care transport,
Adapted for operation aboard aircraft and helicopters,
The device features:
An expansion module that allows the cardiac monitor to be expanded with additional features,
The ability to connect to a defibrillator via the DEFIB SYNCHRO jack,
FREEZE/UNFREEZE function,
Audible and visual alarms with muting capability,
Handle for moving and hanging,
Built-in printer for 60 mm wide thermal paper,
7" liquid crystal monitor,
English language menu,
Keypad on the control panel in the following language: English,
The expansion module has SNAP SHOT and PRINT TREND options,
The device has traces of use in the form of a scratch on the screen and small cracks and abrasions of the case,
Equipped with the following modules with full wiring and sensors:
ECG / ECG RESP - 3 - lead wire,
SOP
2
- New silicone sensor,
NIBP - New cuff with air tubing,
TEMP 1/2 - No sensor included
(can be retrofitted at additional cost),
INV. BP P1/P2 - Missing measurement line
(can be retrofitted at additional cost),
CO
2
- No measuring line
(can be retrofitted at additional cost),
Included:
English manual (PDF),
19 V power supply,
Dimensions: 20.5 x 19.5 x 25 cm,
Weight of the set: 7.5 kg,
Cat 03,
Has a current Technical Passport issued,
Warranty:
6 months for domestic market - Poland,
3 months for the international market,
Possibility of financing (Poland only): Installments, Leasing, Loan,



