



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MIDMARK RITTER 204 gynecological chair (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Made in the U.S.A. (2012),
Medically certified blue upholstery (upholstery can be replaced with new one from our color chart for an additional fee),
The Ritter 204 chair
provides comfort for the patient, ensuring a safe and secure experience during their visit. Drawers on three sides allow easy access to
necessary tools and accessories,
The modern design perfectly complements contemporary furnishing trends in medical, gynecological, and aesthetic medicine offices.
The chair's drawers are made of smooth plastic with rounded edges, easy to keep clean, and consist of a single piece. They ensure a tight seal
against spills inside the drawer.
The enlarged footrest allows the patient to feel safe and secure when getting in and out of the chair, requiring minimal effort.
Features:
Removable seat upholstery,
One-piece drawers,
Large footrest,
Spacious side drawers,
Easy access + capacity,
Extendable stirrups,
Side patient handles,
Dimensions:
Backrest: 94 x 54 cm,
Seat: 70 x 54 cm,
Total dimensions: 85 x 148 cm,
Total length: 148 cm,
Maximum load capacity: 225 kg,
Weight: 110 kg,
Currently inspected and ready for use.
Issued Technical Passport (Service Report) valid for 12-month warranty
Warranty:
6 months for the domestic market (Poland),
3 months for the international market
Possibility of extending the warranty for an additional fee to 12 months for the domestic and international markets
Financing options (Poland only): Installments, Leasing, Loan

