



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



SCHMITZ ARCO MEDI-MATIC 114.500 gynecological chair (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Made in Germany,
Power supply: 230 V,
Frequency: 50/60 Hz,
New upholstery with a medical certificate in gray (upholstery can be changed from our color chart for an additional fee),
The ARCO gynecological chair
is designed for gynecological examinations and procedures.
Technical specifications:
Electric seat tilt adjustment within a range of 50° using a foot controller,
Freely adjustable leg rest abduction within a wide range,
Sliding waste bowl,
Includes:
Power cord,
Wired foot controller,
Sliding waste bowl dimensions: 29 x 32 cm,
Dimensions:
Headrest diameter: ø30 cm,
2-segment backrest:
Upper: 56 x 30 cm,
Lower: 57 x 55 cm,
Maximum load capacity: 180 kg,
Weight: 115 kg,
Updated inspection and ready for operation,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6 months for the domestic market (Poland),
3 months for the international market,
Possibility to extend the warranty for an additional fee to 12 months for the domestic and international markets,
Financing options (Poland only): Installments, Leasing, Loan,



