



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **INVACARE PERFECTO 2 Oxygen Concentrator (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Actual product photos,  
Manufactured by P.R.C.,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 280 W,  
The Perfecto 2 oxygen concentrator  
continuously concentrates oxygen by separating the oxygen contained in atmospheric air from the nitrogen. The operation of the oxygen  
concentrator has no impact on the oxygen content in the air in the room where the device is installed.  
Advantages:  
Very easy to use - ideal for the elderly,  
Quiet and compact - suitable even for demanding users,  
Suitable for continuous operation,  
Mobile on wheels,  
Technical parameters:  
Oxygen concentration level: 94.5%,  
Flow range: 0.5 - 5 l/min (steplessly adjustable using a knob),  
Maximum outlet pressure: 5 psi ( $\pm$  0.5 psi),  
Oxygen sensor (flowmeter) equipped with alarms warning of lower than recommended oxygen concentration,  
Equipped with a mileage counter,  
Noise level: 40 dBA,  
Includes:  
Nasal cannula or oxygen mask with 2.1 m tubing - customer's choice,  
Humidifier (bottle),  
User's manual in English and French, IT, GER, PL, CS,  
Power cord,  
Dimensions: 38 x 30 x 58 cm,  
Weight: 22 kg,  
Updated inspection and ready to use,  
Issued Technical Passport (Service Report) valid for 12 months,  
Warranty:  
6 months for the domestic market (Poland),  
3 months for the international market,  
Possibility to extend the warranty for an additional fee to 12 months for the domestic and international markets,  
Financing options (Poland only): Installments, Leasing, Loan,

