



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Hospital bed (ICU) FAMED ŻYWIEC LE-12 (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- Polish production,
- Power supply: 230 V,
- Emergency power supply: via 4 new built-in batteries 6 V / 1.2 Ah,
- Frequency: 50/60 Hz,
- Power: 350 VA,
- **Technical specifications:**
  - **Electric adjustment via handheld remote control:**
    - Heights: from 45 to 85 cm,
    - Back section elevation angle: 0 / 74°,
    - Thigh segment elevation angle: 0 / 30°,
    - Lower leg segment elevation angle: 24°,
  - **Electric adjustment via central remote control:**
    - Heights: from 45 to 85 cm,
    - Back section elevation angle: 0 / 74°,
    - Thigh segment elevation angle: 0 / 28°,
    - Back section elevation angle 74° and lower leg: 28°,
    - Trendelenburg position: 15°,
    - Anti-Trendelenburg Position: 17°,
    - Sitting position,
    - Function that sets the highest height position,
    - Zeroing function to the zero position,
    - Possibility to block individual functions using switches on the remote control,
  - Built-in CPR safety system,
  - Manually retractable side rails on both sides,
  - Rolling bed with built-in central brake and guidance system,
- **Included in the set:**
  - Power cord,
  - Wired handheld remote control,
  - Central wired remote control,
- Dimensions:
  - Total: 96.5 x 228 x 45/85 cm,
- Maximum load capacity: 250 kg,
- Weight: 175 kg,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
  - **6 months for the domestic market (Poland),**
  - **3 months to the international market,**
  - **Possibility of extending the warranty for an additional fee up to 12 months for the domestic and international markets,**
- Financing options (Poland only): Installments, Leasing, Loan,

