



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



EDAN F9 CTG device (Ex-Display)

New (ex-display),
Actual product photos,
Produced by P.R.C., (2022),
Power supply: 230 V,
Backup power: via a new built-in 14.8 V / 5 Ah battery,
Frequency: 50/60 Hz,
The EDAN F9 ECG
is a solution for all areas of fetal monitoring, used before, during, and after labor. Integrated monitoring of twin pregnancy heart rate, uterine contractions, and fetal movements (manually via a marker or automatically) makes the F9 monitor the most cost-effective and flexible solution for all fetal monitoring needs, without compromising quality, precision, or ergonomics.

Specifications:

- Computer-assisted CTG analysis,
- Twin pregnancy monitoring,
- 24-hour cardiotocography recording,
- 150/152-mm paper for printing test results,
- 1 / 2 / 3 cm/min. - Real print speed,
- Fast printout for saved tests,
- Battery ensuring 2-hour CTG operation without mains power,
- 12.1-inch high-resolution LCD touchscreen with the ability to raise and lower, allowing information such as patient demographics to be entered using the on-screen QWERTY keyboard,
- Prevents monitoring of one fetus with two transducers in the case of twins,
- Water-resistant transducers,
- Compatible with fetal telemetry,
- Stand-alone operation or use in an obstetric monitoring system,
- FHR ultrasonic wave intensity (mW/cm²): 3.6,
- FHR range (bpm): 50 - 240,
- Automatic fetal movement detection,
- Multiple pregnancy monitoring,

The CTG device is equipped with:

- MT-803 EDAN medical cart on a mobile base with 4 braking casters,
- Built-in handle for easy movement trolley,
- Sliding upper and lower drawers for accessories,
- Language selection: ENG, CHIN, ESP, PL, GER, TUR, IT, RUS, GREK, HRV, PT, CS, SWE, DUN, FI, NOR, DEN, BZ,

Includes:

- EDAN F9 MS3-109301 Ultrasound Transducer (B),
- EDAN F9 MS3-109301 Ultrasound Transducer (C),
- EDAN F9 MS3-31527 Ultrasound Transducer (B),
- 3 x transducer straps,
- CTG gel,
- Instruction manual in Polish,
- Power cord,

Dimensions:

- EDAN F9 CTG: 34 x 33 x 12 cm,
- EDAN MT-803 trolley:
- Top shelf: 44.5 x 45 cm,
- Bottom shelf: 37 x 31 cm,
- Sliding top shelf: 22 x 25 cm,
- Sliding bottom shelf: 29 x 31 cm,
- Total dimensions: 44.5 x 50 x 92.5 cm,
- Total weight: 30 kg,

Valid inspection and ready for operation,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty: 6 months for the domestic market (Poland),
3 months for the international market,
Possibility to extend the warranty for an additional fee to 12 months for the domestic and international markets,
Financing options (Poland only): Installments, Leasing, Loan,

