



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



PHILIPS INTELLIVUE MP30 ANESTHESIA patient monitor/cardiometer (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Made in Germany (2012),
Power supply: 230 V,
Battery backup: via 2 built-in batteries – batteries not included
(new batteries can be retrofitted for an additional fee)

Frequency: 50/60 Hz,

The IntelliVue MP30

is a patient monitor that monitors multiple vital signs.

A 10-inch TFT LCD panel display with SVGA resolution.

Operation via a touchscreen and a navigation joystick allows direct access to many functions using simple commands.

A comfortable handle and durable housing make it easy to carry.

Measurement trends can be presented graphically, combining up to three curves on a single graph, allowing for easy identification of changes in the patient's physiological condition.

Parameter trends Rapid-variable data can be updated per beat, and up to four high-resolution trend segments can be displayed.

IntelliVue Technology - access to information from various systems in the hospital network - pharmacy, laboratory, ECG management, radiology PACS, server-based hospital systems, and others - and deliver it to the patient monitor.

No separate hard drive or fan - IntelliVue monitors are small, quiet, stable, and start quickly.

The navigation knob makes it easy to enter information and navigate the screen.

Acoustic and visual alarms.

LAN.

Patient monitor equipped with the PHILIPS IntelliVue X2 measurement module:

ECG/RESP - New 3-lead cable.

SPO2 - New silicone sensor.

NBP - New cuff with air tube.

TEMP - Temperature connector socket.

(sensor not included, can be retrofitted for an additional fee.)

Menu in Language: Polish

Includes:

Instruction manual in English (PDF),

Power cord,

New 3-lead ECG cable,

New silicone SPO2 sensor,

New NIBP cuff with air tube,

Dimensions: 21 x 34 x 28 cm,

Weight: 6.5 kg,

Cat 04,

Updated inspection and ready to use,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6 months for the domestic market (Poland),

3 months for the international market,

Possibility of extending the warranty for an additional fee to 12 months for the domestic and international market,

Financing options (Poland only): Installments, Leasing, Loan,

