



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



SMITHS LEVEL 1 EQUATOR EQ-5000 Patient Warmer (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- Made in the USA,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- **The EQUATOR** convection heating system provides controlled convection heating in a compact package,
- Electronically controlled temperature regulation using a thermoelectric sensor,
- It can prevent all clinical effects of hypothermia, improve postoperative comfort of patients and shorten the stay in the recovery room,
- The device can be free-standing or mounted on an IV stand, bed frame or optional wheeled trolley,
- **Effective temperatures:**
 - High: 44° C +/- 1° C,
 - Average: 40° C +/- 1° C,
 - Low: 36° C +/- 1° C,
- Time to change the average temperature from 20 ° C to 36 ° C: 7 minutes,
- **Equipment placed on a dedicated mobile cart,**
- **Included in the set:**
 - Power cord,
 - New biological filter,
 - Pipe with hot air supply,
- Dimensions: 24 x 30 x 26 cm,
- Weight: 10 KG,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months to the international market,**
 - **Possibility of extending the warranty for an additional fee up to 12 months for the domestic and international markets,**
- Financing options (Poland only): Installments, Leasing, Loan,

