



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



#### **ALVO MEDICAL ETIUDA 4.02 operating table (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Illustrative photos,  
Made in Poland,  
6-segment,  
Power supply: 230 V, and 2 x 12V 7A batteries (new batteries),  
Frequency: 50-60 Hz,  
Power: 350 VA,  
Electrical adjustment of settings using a wired remote control:  
Head orientation selection,  
Table "0" function,  
Up-down range: 75-108 cm,  
Vertical tilt - Trendelenburg, Reverse Trendelenburg position, +35°/-35°,  
Lateral tilt: +20°/-20°,  
Manual adjustment of settings:  
Manually adjustable foot segment -100°/+20°, via gas springs, divided and laterally retracted 90°,  
Manually adjustable headrest +35°/-20° via gas springs,  
Tabletop slide range: 30 cm,  
Manually adjustable backrest +80°/-20° via gas springs,  
Manual lumbar adjustment range: 11 cm (crank),  
Manual brake,  
Features:  
X-RAY transparency,  
Rails for mounting additional accessories,  
Bed dimensions:  
Width: 50 cm,  
Length: 210 cm,  
Thickness: 6 cm,  
Base made of stainless steel,  
Base dimensions: 58 x 99 cm,  
Overall dimensions: 58 x 210 x 75/155 cm,  
Permissible load: 250 kg,  
Mobile on casters with a central brake,  
Currently inspected and ready for use ready for operation,  
Issued Technical Passport (Service Report) valid for 12 months,  
Warranty:  
6 months for the domestic market (Poland),  
3 months for the international market,  
Possibility of extending the warranty for an additional fee to 12 months for the domestic and international markets,  
Financing options (Poland only): Installments, Leasing, Loan,  
DESCRIPTION AND ILLUSTRATIVE PHOTOS OF THE TABLE

