



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ADEC 500 dental unit (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Made in the U.S.A.,
Power supply: 230 V,
Frequency: 50/60 Hz,
Medically certified orange upholstery (upholstery can be replaced with a new one from our color chart for an additional fee),
The ADEC 500 dental unit
is a compact, ergonomic dental unit with an electromechanical drive,
The operator's table is located on a pantograph arm:
5 upper silicone sleeves:
Air-water syringe (three-way),
New polymerization light
COXO D3-686-1B+
New scaler with ENDO function and LED light
WOODPECKER UDS-N3LED
(autoclavable at 135°),
2x new sleeves Silicone
MINDWEST
with light -
A turbine and micromotor (with light) can be retrofitted for an additional fee.
Smooth control of pneumatic tips via foot control.
Freely extendable, manually operated operator tray.
Halogen operator lamp
ADEC 17 V / 95 W
with free vertical and horizontal adjustment -
Lamp can be replaced with an LED lamp for an additional fee.
Electric seat adjustment via foot control:
Seat height adjustment: 48 - 75 cm.
Backrest tilt: -15° / 0° / 70° (automatic seat backrest synchronization).
4 patient functions programmed on the foot control.
The chair can be reclined to the Trendelenburg position.
Freely manually adjustable headrest height and tilt.
A new wet pump can be retrofitted. Suction pump and compressor available for an additional fee.
Includes:
Air-water syringe (three-way),
New polymerization lamp
COXO D3-686-1B+,
New scaler with ENDO function and LED light
WOODPECKER UDS-N3LED
(autoclavable at 135°),
2x new silicone sleeves
MINDWEST
with light,
Dimensions:
Headrest: 20 x 18 cm,
Backrest: 65 x 45 cm,
Seat: 58 x 120 cm,
Total dimensions: 90 x 130 x 194 cm,
Weight: 170 kg,
Updated inspection and ready for use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty: 6 months for the domestic market (Poland),
3 months for the international market,
Possibility of extending the warranty for an additional fee to 12 months for the domestic and international markets,
The warranty period does not cover light bulbs and other light sources installed in the device,
Financing options (Poland only): Installments, Leasing, Loan,

