



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **MIDMARK 630-002 gynecological treatment chair (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Actual product photos,  
Made in the U.S.A.,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Medically certified gray upholstery – the upholstery shows signs of wear (upholstery can be replaced with a new one from our color chart for an additional fee),  
MIDMARK 630-002  
is a treatment chair designed to make patients feel comfortable and safe.  
Thanks to its low height and multi-position adjustments, the chair facilitates easier patient positioning at the touch of a button, reducing the likelihood of fatigue, back and neck strain.  
Smooth electric adjustment via a wired foot switch and handheld remote control:  
Footrest angle from 0° to 90°,  
Seat angle from 0° to 30°,  
Backrest tilt angle from 0° to 85°,  
Seat height range: 62 - 105 cm,  
Trendelenburg position adjustable within 30°,  
4 programmable settings,  
Manual headrest height and tilt adjustment,  
Features:  
Manually adjustable headrest,  
Extending, folding, and side-adjustable stirrups,  
Extending treatment bowl,  
Paper towel compartment with a width of 45 cm,  
Side rails for mounting additional accessories,  
Includes:  
Wired foot control,  
Wired handheld remote control,  
English user manual (PDF),  
Power cord,  
Dimensions:  
Headrest: 30 x 32 cm,  
Backrest: 68 x 58 x 10 cm,  
Seat: 68 x 52 cm x 10 cm,  
Footrest: 42 x 25,  
Base: 52 x 100 cm,  
Maximum stirrup spacing: 114 cm,  
Overall chair length in the lying position: 185 cm,  
Overall dimensions: 116 x 74 x 145 cm,  
Maximum load capacity: 204 kg,  
Weight: 195 kg,  
Updated inspection and ready for use,  
Issued Technical Passport (Service Report) valid for 12 months,  
Warranty: 6 months for the domestic market (Poland),  
3 months for the international market,  
Possibility of extending the warranty for an additional fee to 12 months for the domestic and international market,  
Financing options (Poland only): Installments, Leasing, Loan,

