



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Ophthalmic ultrasound device ECHOSON DESMIN H-USO (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Made in Poland,
Power supply: 230 V,
Frequency: 50/60 Hz,
Specifications:
Built-in monitor,
Display:
Grayscale: 256 levels,
Image matrix: 512 x 512 px.,
LCD: 10.4",
Eye tests:
Lentilated,
Acrylic lens,
PMMA lens,
Cataract,
Lentless,
Other lens
Operation algorithm: fast / manual / advanced,
Amplitude level (%): 60 / 70 / 80 / 90 / 100,
Menu language: English,
Includes:
Power cord,
ECHOSON OA 12 MHz probe,
Probe ECHOSON OB 12 MHz,
MITSUBISHI P93E Video Printer,
Overall dimensions: 46 x 51 x 100 cm,
Weight: 8.2 kg,
Updated inspection and ready to use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6 months for the domestic market (Poland),
3 months for the international market,
Possibility to extend the warranty for an additional fee to 12 months for the domestic and international markets,
Financing options (Poland only): Installments, Leasing, Loan,

