



DEVILBISS HEALTHCARE 525KS Oxygen Concentrator (Ex-Display)

New (ex-display),
Actual product photos,
Made in the U.S.A.,
Power supply: 230 V,
Frequency: 50/60 Hz,

The Devilbiss Compact 525KS oxygen concentrator is a reliable device designed for oxygen therapy in various settings, such as homes, nursing homes, and assisted living facilities.

It is especially recommended for patients suffering from COPD, cardiovascular disease, and lung diseases.

Features:

Very quiet operation: The noise level is only 40 dBA, ensuring comfortable use without disruptions.

Easy operation: The intuitive control panel allows for easy use.

High oxygen concentration: 94% at a flow rate of 0.5 to 5 liters per minute ensures effective therapy.

Adjustable flow: The oxygen flow can be set from 0 to 5 liters per minute, tailored to the patient's individual needs.

Alarm system: A visual and audible low output alarm ensures safety. Use,

Mobility: Equipped with wheels, allowing for easy movement of the device.

Includes:

Nasal cannula or oxygen mask with a 2.1 m tube - to be chosen by the customer.

Condensate container.

Power cord.

Dimensions: 34 x 30 x 62 cm.

Weight: 16.4 kg.

Updated inspection and ready for operation.

Issued Technical Passport (Service Report) valid for 12 months.

Warranty:

6 months for the domestic market (Poland),

3 months for the international market.

Possibility of extending the warranty for an additional fee to 12 months for the domestic and international markets.

Financing options (Poland only): Installments, Leasing, Loan.

