



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Polometer (perimeter) CARL ZEISS HUMPHREY SYSTEM 720i (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- Made in the USA,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Bulb: Halogen 12V/30W,
- **Technical specifications:**
 - Field of view: 90°,
 - Maximum stimulus intensity: 10,000 ASB,
 - Stimulus duration: 200 ms,
 - Stimulus wavelength: Broadband visible light,
 - Canopy illumination: 31.5 ASB,
 - Test distance: 30 cm,
 - **Threshold Test Library:**
 - 24-2, 30-2, 10-2, Spot,
 - 60-4,
- **Threshold Testing Strategies:**
 - SITA Standard, SITA Fast, Full Threshold, FastPac,
- **Screening Test Library:**
 - Central Field Test Patterns: C40, C76, C80,
 - Peripheral field test patterns: P-60, FF-81, FF-120, FF-135,
 - Screening Test Strategies: Two Zone, Three Zone, Quantify Defects,
 - Screening Test Modes: Two Zone/Three Zone and Defect Determination, Age-Related Level and Threshold Modes, Single Intensity,
- **Special test library:**
 - Esterman Monocular,
 - Esterman's binoculars,
 - Superior Room 36,
 - Superior Room 64,
 - Stimulus / Background Color: white on white,
 - General test features: Goldmann stimulus sizes: III,
- User-defined test storage,
- It has a built-in 11" touchscreen monitor.
- Electric chin rest adjustment: up-down and sideways,
- **Possibility to select menu in: EN, GER, ES, FR, IT, JP, PT, SE,**
- **The whole is placed on a mobile ophthalmological table MDT MD-22 with the possibility of electric adjustment in the range of: 75 - 95 cm,**
- **Included in the set:**
 - User manual in English (PDF),
 - Power cord,
 - Patient button,
 - D-SUB25 printer cable,
 - HD-LASERJET M60 printer,
- Overall dimensions: 86 x 43 x 125 / 145 cm,
- Weight: 85 KG,
- Cat 04,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for a period of 12 months,**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months to the international market,**
 - **Possibility of extending the warranty for an additional fee up to 12 months for the domestic and international markets,**
 - **The warranty period does not cover light bulbs and other light sources installed in the device,**
- Financing options (Poland only): Installments, Leasing, Loan,

