



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ERBE ERBOKRYO AE Cryosurgery Apparatus (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- German production,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 40 W,
- **Application:**
 - Ophthalmology,
 - Dermatology,
 - Otolaryngology,
 - Proctology,
 - Veterinary surgery,
 - Plastic surgery,
 - Gynecology,
 - Thoracic surgery,
- **Specification:**
 - Temperature setting range: From -20 °C to -70 °C,
 - Defrosting temperature: From 8 °C to 15 °C,
 - **Cryogenic gas:**
 - Gas type: Nitrous oxide (N₂O) or carbon dioxide (CO₂),
 - Working pressure: 40 – 60 bar,
 - Gas consumption "freezing": 36 g/min.,
 - Amount of gas produced: 40 – 60 l/min.
- **Included in the set:**
 - ERBE probe 20716-087 Ø 3.0 mm.,
 - Hose to the supply bottle,
 - Foot manipulator,
 - Power cord,
- Overall dimensions: 47 x 35 x 10 cm,
- Weight: 10 KG,
- Cat 03,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months to the international market,**
 - **Possibility of extending the warranty for an additional fee up to 12 months for the domestic and international markets,**
- Financing options (Poland only): Installments, Leasing, Loan,

