



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



CARL ZEISS HUMPHREY FIELD ANALYZER 745i (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Made in the U.S.A.,
Frequency: 50/60 Hz,
Halogen bulb: 12 V / 20 W,
The Humphrey HFA 745i
is an automated circuit designed to identify visual field defects for screening, monitoring, and assisting in the diagnosis and treatment of eye diseases such as glaucoma and related neurological disorders.
Ergonomic design ensures maximum comfort, access, and versatility.
Equipped with an automatic head and eye tracking and alignment system.
Test types:
Static,
Kinetic,
Blue on Yellow,
Specifications:
Field of view: 89°
Maximum stimulus intensity: 10,000 ASB,
Stimulus duration: 200 ms,
Wavelength Stimulus Distance: 30 cm,
Cup Illumination: 31.5 ASB,
Testing Distance: 30 cm,
Features:
Threshold Test Library:
24-2, 30-2, 10-2, Spot,
60-4, Nasal Step,
Threshold Test Strategies:
SITA Standard, SITA Fast, Full Threshold, FastPac,
SITA-SWAP,
Screening Test Library:
C40, C64, C76, C80, C-Armal,
Peripheral Test Patterns,
Screening Test Modes:
Age Correction,
Threshold-Related, Single Intensity,
Special Test Library:
Social Security Disability, Monocular, Binocular,
Superior 36, 64,
Custom Testing,
Fixation Control:
Heijl Blind Spot Monitor Krakau,
Video Eye Monitor,
Gaze Tracking,
Remote Video Eye Monitor Capability,
Stimulus:
White on White,
Red or Blue on White,
Blue on Yellow (SWAP),
General Test Features:
Goldmann I-V,
Lower Threshold Test,
Software Features:
STATPAC 2 - Single Field Analysis,
Glaucoma Hemifield Test (GHT),
Visual Field Index (VFI),
Guided Progression Analysis (GPA),
Serial Field Review,
Network,
FORUM Connectivity,
DICOM Connectivity,
EasyConnect RCT / HFA-NET Pro. Features:
A built-in touchscreen and menu-driven interface simplify operation,
Electric chin rest adjustment: up-down and sideways,
Built-in floppy disk drive,
Menu language selection: ENG, GER, ESP, FR, IT, JP, PT, SWE,
Includes:
English user manual (PDF),
Power cord,
LEXMARK laser printer,
Wired mouse,
Wired keyboard,
Patient button,
D-SUB25 printer cable,
Dimensions:
Perimeter: 55 x 55 x 60 cm,
Printer: 40 x 40 x 26 cm,
Weight: 60 kg,
Updated inspection and ready for use,
Issued Technical Passport (Service Report) valid for 12-month warranty
Warranty:
6 months for the domestic market (Poland),
3 months for the international market
Possibility of extending the warranty to 12 months for the domestic and international markets for an additional fee
Warranty period does not cover light bulbs and other light sources installed in the device
Financing options (Poland only): Installments, Leasing, Loan

