



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ERBE medical equipment trolley (Reconditioned)

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Actual product photos,

Made in Germany,

ERBE mobile medical cart designed for medical equipment and cryosurgery procedures,

Features:

Stable and durable construction,

Tabletop coating adapted for easy cleaning,

Equipped with a rear mounting for cylinders,

Sliding drawer for equipment,

Assistant placed on 4 casters (2 with built-in brakes),

Dimensions:

Top shelf: 56 x 38 cm,

Sliding drawer: 48 x 32 x 8 cm,

Total dimensions: 64 x 46 x 82 cm,

Weight: 15 kg,

Updated inspection and ready for operation,

Issued Technical Passport (Service Report) valid for 12-month warranty

Warranty:

6 months for the domestic market (Poland),

3 months for the international market

Possibility of extending the warranty for an additional fee to 12 months for the domestic and international markets

Financing options (Poland only): Installments, Leasing, Loan

