



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



SCHMITZ MEDI-MATIC gynecological chair 115.755 (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- German production,
- Frequency: 50/60 Hz,
- **The SCHMITZ MEDI-MATIC gynecological chair is used to perform gynecological and urological examinations and procedures,**
- **Medically certified upholstery in orange (the upholstery can be replaced with a new one from our color chart for an additional fee),**
- **Technical specifications:**
 - **Smooth electric adjustment via foot controller:**
 - Heights in the range: 60 – 90 cm,
 - Backrest adjustment angle in the range: $\pm 52^\circ$,
 - Seat up/down angle: $\pm 5^\circ$,
 - Trendelenburg position: 19° ,
 - **Smooth electric adjustment via handheld remote control:**
 - Heights in the range: 60 – 90 cm,
 - Backrest adjustment angle in the range: $\pm 52^\circ$,
 - Seat up/down angle: $\pm 5^\circ$,
 - Trendelenburg position: 19° ,
 - Double-sided footrests in the range of: $\pm 45^\circ$,
- Quick change of the seat position using the foot controller,
- The low starting position of the seat allows for easy and safe entry without the need to use a step,
- Paper roll storage,
- Paper roll hanger, 48 cm wide,
- Pull-out waste bowl from under the seat section,
- Gentle start and ergonomic movement when changing position,
- **Included in the set:**
 - Power cord,
 - Wired foot controller,
 - Wired remote control,
 - A bowl for waste,
 - Headrest,
- **Dimensions:**
 - Bases: 105 x 58 cm,
 - Waste bowls: 31 x 26 cm,
 - Overall: 150 x 135 x 86 cm,
- Maximum load capacity: 200 KG,
- Weight: 165 KG,
- Cat 04,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months to the international market,**
 - **Possibility of extending the warranty for an additional fee up to 12 months for the domestic and international markets,**
- Financing options (Poland only): Installments, Leasing, Loan,

