



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **PENLON SIGMA DELTA (Isoflurane) Inhalant Vaporizer (Reconditioned)**

Reconditioned (used),

Technical condition: very good,

Visual condition: good,

Actual product photos,

Made in the U.K.,

Power supply: 230 V,

Frequency: 50/60 Hz,

**PENLON SIGMA DELTA**

is a safe and precise vaporizer for inhalants.

Specifications:

Flow range: 0.2-15 l/min.,

Concentration adjustment: 0.2%-5%,

Tank capacity: 250 ml,

Compensation: automatic flow and temperature correction.

Control panel: large, easy-to-read dial + anesthetic level indicator.

Infusion type: Key-Fill.

This vaporizer uses the DRAGER DW-2000 mounting system.

Includes:

Instruction manual in English (PDF).

Dimensions: 10 x 19 x 24 cm. cm,

Weight: 4.8 kg,

It has a current inspection and is ready for operation,

A Technical Passport (Service Report) is issued, valid for 12 months,

Warranty:

6 months for the domestic market (Poland),

3 months for the international market,

Possibility of extending the warranty for an additional fee to 12 months for the domestic and international markets,

Financing options (Poland only): Installments, Leasing, Loan,

