



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Videoprinter MITSUBISHI P93E (Reconditioned)

Reconditioned (used),
Technical condition: very good,
Visual condition: very good,
Actual product photos,
Made in Japan,
Power supply: 230 V,
Frequency: 50/60 Hz,
The Mitsubishi P93E
is a fast, high-resolution monochrome thermal printer designed for printing images from video devices, including: from medical diagnostic devices, CCTV cameras, and VCRs with full photographic quality.
Specifications:
Long-lasting print heads,
Thermal printing - true photorealism - 256 grayscale levels,
Roll Type technology - prints on 110 mm wide paper rolls,
Print format: 100 x 75 mm (Normal) and 93.7 x 125 mm (Postcard), as well as 50% Normal and 50% Postcard,
Print time: 3.5 seconds,
Print resolution: 325 dpi,
Image memory for 10 individually selected frames,
Inputs: Composite PAL, signal system: PAL / NTSC video,
Printing in 2 or 4 parts, Multiimage,
Function menu displayed on the LED display,
Panel menu in English,
Printing of setting parameters under the photo,
Power saving mode Paper capacity (260 Normal prints per roll),
Function lock,
Includes:
User manual in English (PDF),
Power cord,
Wired remote printing button,
BNC cable,
Dimensions: 15.5 x 27 x 9 cm,
Weight: 3 kg,
Updated inspection and ready to use,
Issued Technical Passport (Service Report) valid for 12 months,
Warranty:
6 months for the domestic market (Poland),
3 months for the international market,
Possibility to extend the warranty for an additional fee to 12 months for the domestic and international markets,
Financing options (Poland only): Installments, Leasing, Loan,

