



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



### **Surgical diathermy ERBE ICC 80 (INT) (Reconditioned)**

Reconditioned (used),  
Technical condition: very good,  
Visual condition: very good,  
Real product photos,  
German production,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 80 W,  
Smooth power regulation,  
Digital displays of set power,  
Foiled, easy-to-use control panel,  
Has the following features:  
Two types of monopolar cutting - max. 80 W:  
Monopolar cutting,  
CUT CONTROL - optimal, repeatable cutting,  
Three types of coagulation - max. 50 W:  
SOFT COAG - precise control of the coagulation process and prevention of carbonization through coagulation mode,  
FORCED COAG - when fast and efficient coagulation is needed,  
BIPOLAR COAG - bipolar coagulation,  
Included:  
Brand-new monopolar electrode lead with ERGO PEN BOWA handle - two buttons, 4 mm shaft (no electrodes included, possibility of upgrading for an additional fee),  
Brand-new lead with BOWA clip for disposable neutral electrodes (no electrodes included, possibility of upgrading for an additional fee),  
Brand-new bipolar lead, BOWA forceps  
(possibility to upgrade with a bipolar forceps for an additional fee),  
Corded scalpel attachment,  
User manual in English (PDF),  
Power cord,  
Dimensions: 27 x 11 x 27 cm,  
Weight: 2.8 KG,  
Has current inspection and is ready for use,  
Technical Passport (Service Report) issued and valid for 12 months,  
Warranty:  
6 months for the domestic market (Poland),  
3 months for the international market,  
Possibility to extend warranty for an additional fee up to 12 months for domestic and international markets,  
Financing available (Poland only): Installments, Leasing, Loan,

