



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



PHILIPS INVIVO EXPRESSION MRI Patient Monitor/Cardiomonitor (Reconditioned)

- Reconditioned (used),
 - Technical condition: very good,
 - Visual condition: very good,
 - **Real product photos,**
 - Made in the USA (2014)
 - Power supply: 230 V,
 - Emergency power supply: via 4 new batteries 14.8 V / 75 Wh,
 - Frequency: 50/60 Hz,
 - **System equipped with:**
 - **INVIVO EXPRESSION patient monitor on MR cart with built-in anesthetic gas measurement module:**
 - **Invivo Expression** is an MRI patient monitor designed specifically for environments with very strong magnetic fields,
 - The monitor is equipped with wireless connectivity for ECG and SpO2 monitoring,
 - 360-degree alarm light notifies you of an alarm condition when traditional alarms are violated,
 - **Specification:**
 - 12" color LCD screen,
 - Innovative user interface,
 - Resolution 800×600 pixels,
 - Quickly connect accessories,
 - **Monitor integrated with the mobile MR cart:**
 - **Measurement:**
 - 2x invasive blood pressure measurement,
 - 2x NIBP measurement,
 - Low-flow CO2 / N2O / O2 measurement (**built-in module for measuring anesthetic gases**) ,
 - Wireless SPO2 measurement,
 - Wireless ECG measurement,
 - **PHILIPS EXPRESSION MR IP5 monitor:**
 - Reconditioned (used),
 - Technical condition: very good,
 - Visual condition: very good,
 - **Real product photos,**
 - Made in the USA (2020)
 - Power supply: 230 V,
 - Frequency: 50/60 Hz,
 - **The Expression IP5 system** is a solution that enables the display of vital signs through remote viewing and monitoring, without the need for extensive and costly wireless networks to connect the Expression IP5 system to Expression series MR patient monitoring solutions.
 - **Characteristics:**
 - Integral color LCD screen,
 - Intuitive graphical interface,
 - Colored waveform and large numbers,
 - Up to 8-hour battery life with built-in charger,
 - Simultaneous display of parameter, waveform and related values,
 - Visual and audible alarm signals, heart rate signal, system status information,
 - **Specification:**
 - Color VGA liquid crystal display, 640×480 pixels, five-wire touch screen,
 - Diagonal 5.7" inches (14.5 cm),
 - Swipe speed: 25mm/s giving 2.7 seconds of display time,
 - Waveform display mode: fixed record/moving,
 - **Application:**
 - Sedation,
 - Anesthesiology,
 - Pediatrics,
 - Obstetrics,
 - ECG gating,
 - Interventional procedures,
 - Intraoperative procedures,
 - Patient transport in the MR environment,
- **Possibility to select menu in: ENG, ESP, FR, PT, IT, DUN, SWE, NED,**
- **Included in the set:**
 - PHILIPS MR USB Thermal Printer,
 - USB printer cable,
 - Wireless SPO2 measurement module,
 - Wireless ECG measurement module (4-lead),
 - New NIBP cuff with air tube,
 - Battery charging station for patient measurement modules (4 batteries included),
 - PHILIPS IP5 EXPRESSION monitor power supply,
 - Power supply for the INVIVO EXPRESSION measurement kit,
 - Power cord,
- Dimensions:
 - PHILIPS IP5 EXPRESSION monitor: 50 x 30 50 cm
 - INVIVO EXPRESSION monitoring system: 44 x 50 140 cm,
- Total weight: 80 KG,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months to the international market,**
 - **Possibility of extending the warranty for an additional fee up to 12 months for the domestic and international markets,**
- Financing options (Poland only): Installments, Leasing, Loan,

