



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ZOLL R SERIES PLUS Defibrillator (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- Made in the USA,
- Power supply: 230 V,
- Emergency power supply: via 1 new built-in battery 10.8 V / 5.8 Ah - **(not included) can be retrofitted for an additional fee,**
- Frequency: 50/60 Hz,
- **Menu in: English,**
- **The R Series defibrillator** is widely used in hospitals, especially in cardiology and intensive care units,
- Thanks to the use of advanced ZOLL technologies, it ensures high effectiveness of electrotherapy (defibrillation, cardioversion, pacing) and enables real-time monitoring of CPR quality,
- Provides audio and visual feedback to assist the rescuer in consistently providing high-quality CPR,
- The visual perfusion index (PPI) indicates whether the combination of compression rate and depth is consistent with the guideline recommendations,
- **Characteristics:**
 - **Real CPR Help for ALS** - Provides real-time feedback on compression depth and speed,
 - **OneStep System** - simplifies defibrillation, pacing and monitoring with one set of electrodes,
 - **Code-Readiness Testing System** - ensures that the device is always ready for use in case of emergency,
 - **Rectilinear biphasic defibrillation technology** - provides more effective defibrillation with less energy consumption,
 - **Defibrillation Advisory** - complies with AHA/ERC guidelines for shock delivery,
 - **Built-in test port** - allows for hassle-free system inspection and maintenance,
 - **Compact Flash card slot** - allows for data storage and easy transfer of critical event logs,
- **Technical specifications:**
 - LCD liquid crystal display: 6.5" inch,
 - **Power selection:**
 - Adult: 1 / 2 / 3 / 4 / 5 / 6 / 7 / 8 / 9 / 10 / 15 / 20 / 30 / 50 / 75 / 100 / 120 / 150 / 200 J,
 - Child: 1 / 2 / 3 / 4 / 5 / 6 / 7 / 8 / 9 / 10 / 15 / 20 / 30 / 50 / 70 / 85 / 100 / 120 / 150 / 200 J **(OneStep pediatric electrodes only)** ,
 - **CPR Monitoring:**
 - Compression depth: From 1.9 cm to 7.6 cm (+/- 0.6 cm),
 - 50 to 150 compressions per minute,
 - Heart rate measurement range: From 30 to 300 BPM,
 - Built-in audible and visual alarm system,
 - Built-in printer with paper width: 9 cm (print speed 25 mm/s),
- **Included in the set:**
 - User manual in English (PDF),
 - Power cord,
 - 3-lead ECG cable,
 - SPO2 cable,
 - ZOLL EtCO2 Sensor,
 - Set of hard spoons **(not included - can be retrofitted for an additional fee),**
- Dimensions: 35 x 34 x 28 cm,
- Total weight: 7 KG,
- Cat D153,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months to the international market,**
 - **Possibility of extending the warranty for an additional fee up to 12 months for the domestic and international markets,**
- Financing options (Poland only): Installments, Leasing, Loan,

