



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Roller mixer COLE-PARMER RS-200 (Reconditioned)

Refurbished (used),

Technical condition: very good,

Visual condition: very good,

Real product photos,

Made in U.K.,

Power supply: 230 V,

Frequency: 50/60 Hz,

Power: 50 W,

Roller mixer

COLE-PARMER

provides gentle yet efficient mixing, ideal for mixing blood samples, viscous substances and suspensions liquid-solid where minimal aeration is required,

RS-200 roller mixers are now equipped with detachable side panels. The side panels provide the user with full assurance that the sample will be effectively mixed during continuous operation, and sample bottles will not slide off the mixer during operation,

Features:

Safe for use in coolers down to 4°C (39°F) or in incubators up to 40°C (104°F) and humidity up to 80%,

Continuous and quiet operation,

Side panels ensure sample bottles remain on the shaker,

Specification:

Number of rollers: 6,

Speed: 33 rpm,

Amplitude: 16 mm,

Roller length: 34 cm,

In the box:

Power cord,

Overall dimensions: 59.5 x 25.5 x 11 cm,

Weight: 5.1 KG,

Includes current inspection and is ready for use,

Technical Passport (Service Report) issued for 12 months,

Warranty:

6 months for the domestic market (Poland),

3 months for the international market,

Possibility to extend warranty for an additional fee to 12 months for domestic and international markets,

Financing option (Poland only): Installments, Leasing, Loan,

