



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **LEICA M620 F20 surgical microscope (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real product photos,**
- German production, (2022),
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Main light source: Halogen 12V / 50W,
- Spare light source: Halogen 12V / 50W,
- **The LEICA M620F20** is a high-quality ophthalmic microscope with a large depth of field. Its advantages include comfortable operation, a high-quality light source, and, above all, excellent optics.
- **Electrically controlled microscope functions:**
  - Smooth adjustment of magnification (zoom) and focus (focus) using the foot controller and control panel,
  - Adjusting the focusing and zoom speed,
  - Control unit with LCD display, allowing you to save configuration settings for several users,
- **Specification:**
  - **Total microscope magnifications:**
    - Doctor: from 3x.9 to 23.2x,
    - Assistant: from 3.1x to 18.8x,
  - Doctor's glasses magnification: x12.5
  - Assistant eyepiece magnification: 10x,
  - Doctor's view with adjustable angle: from 5° to 25°,
  - Assist preview with fixed tilt angle: 45°,
  - Lens focal length f=200 mm. APO,
  - Head tilt mechanism: 5° / -5°
  - XY module reset function - return to the initial position,
  - Focus module reset function - return to the initial position,
  - IR protection filter,
  - UV protection filter: KV420,
  - A module with a spare bulb built into the microscope head,
- Included in the set:
  - Wireless 12-function foot controller,
  - VIDEO recorder,
  - Set of handles (3 pcs.)
- Dimensions:
  - Total height:
    - Minimum: 190 cm,
    - Maximum: 222 cm,
  - Lens distance from ground:
    - Minimum: 78 cm,
    - Maximum: 138 cm,
  - Maximum distance from the lens to the column: 137 cm,
  - Base: 61 x 61 cm,
- Weight: 210 KG,
- **It has a current inspection and is ready to work,**
- **Issued Technical Passport (Service Report) valid for 12 months,**
- **Guarantee:**
  - **6 months for the domestic market (Poland),**
  - **3 months to the international market,**
  - **Possibility of extending the warranty for an additional fee up to 12 months for the domestic and international markets,**
  - **The warranty period does not cover light bulbs and other light sources installed in the device,**
- Financing options (Poland only): Installments, Leasing, Loan,

