



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



NDS ENDOVUE LCD BC-WU24-A1417 Medical Monitor (Reconditioned)

Refurbished (used),
Technical condition: very good,
Visual condition: very good,
Real product photos,
Manufactured in USA,
Power: 230 V,
Frequency: 50/60 Hz,
NDS ENDOVUE

is a monitor capable of displaying multiple imaging modes, designed for minimally invasive surgery. The 24-inch surgical monitor NDS ENDOVUE is characterized by a sturdy, lightweight and flexible construction, enabling easy installation, and patented factory color calibration guarantees reliable, consistent and realistic image reproduction,

This is an ideal solution for professionals providing image quality comparable to high-end endoscopic displays,

Technical specifications:

Display: Active Matrix TFT IPS,

Screen: 24" inches,

Resolution: 1920 x 1200 px. (WUXGA),

Luminance: 300 cd/m² (Max), 250 cd/m² (typical at 6500K),

Contrast: 1000:1,

Maximum viewing angle: 178°,

Response time: 14 ms,

Connections:

DVI / HD-SDI / HD-GBS / RGBS / VGA / S-VIDEO / COMPOSITE,

Included:

User manual in English (PDF),

Power supply with power cord,

Dimensions: 60 x 10 x 40 cm,

Weight: 15 KG,

Has current inspection and is ready for use,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6 months for domestic market (Poland),

3 months for international market,

Option to extend warranty for an additional 12 months for domestic and international markets,

Financing option (only Poland): Installments, Leasing, Loan,

