



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Olympus OFP-2 Flushing Pump (Reconditioned)**

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Real product photos,

Manufactured in the U.K., (2014),

Power supply: 230 V,

Frequency: 50/60 Hz,

Power: 100 VA,

Washing pump OFP-2

provides a strong water jet that can be easily controlled. The device can be connected to the working channel and auxiliary water supply channel of endoscopes from Olympus such as: EVIS, ACTERA, EXERA, LUCERA, SMARTAGE and EUS,

Application:

Improving diagnostics and visualization of the treatment field during endoscopic examinations and endotherapy through effective washing of the stomach and colon mucosa,

Facilitating the use of an endoscopic ultrasound probe by quickly filling the organ with water for an endoscopic USG examination,

Specification:

Microprocessor control ensures accurate and precise water dosing,

Ergonomically designed placement of control buttons, ensuring easy operation of the device,

Bright LED display for easy verification of selected settings,

Construction easy to keep clean and hygienic,

Overpressure protection in Olympus endoscope channels,

Standby mode allows quick and safe replacement of the water feeding tube without using the main power button,

The pump head automatically shuts off when opened, preventing accidental injury to the operator,

Automatic timer turning off the pump after 20 seconds of operation allows safe water delivery to the patient's body,

Maximum water flow:

Tool channel: 700 to 750 ml/min,

Auxiliary channel: 230 ml/min,

In the set:

English-language user manual (PDF),

Foot pedal attachment,

OLYMPUS MAJ-1603 container,

Power cord,

Dimensions: 20 x 17 x 39 cm,

Weight: 4 KG,

Has current inspection and is ready for use,

Technical Passport (Service Report) issued, valid for 12 months,

Warranty:

6 months on the domestic market (Poland),

3 months on the international market,

Warranty extension available for an additional fee up to 12 months for domestic and international markets,

Financing options (Poland only): Installments, Leasing, Loan,

