



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



SOEHNLE baby scale (Reconditioned)

Refurbished (used),
Technical condition: very good,
Visual condition: very good,
Real product photos,
Made in Germany,
Power supply: 230 V,
Frequency: 50/60 Hz,
Electronic scale for weighing infants
SOEHNLE
ensures proper and trouble-free operation during weighing newborns,
Technical specification:
Load capacity: every 10 g,
Maximum newborn weight: 20 KG,
Display: LCD,
Working surface dimensions: 52 x 20 cm,
Maximum load: 20 KG,
Weight: 3 KG,
Has current inspection and is ready for operation,
Technical Passport (Service Report) issued for 12 months,
Warranty:
6 months for the domestic market (Poland),
3 months for the international market,
Warranty extension available for an additional charge up to 12 months for domestic and international markets,
Financing options (Poland only): Installments, Leasing, Loan,

