



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Operating lamp FAMED BH-132 (Reconditioned)**

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Real product photos,

Polish production,

Power supply: 230 V,

Frequency: 50/60 Hz,

Bulb: 24 V / 65 W,

Operating lamp FAMED BH-132

is used as a light source in hospitals, clinics, veterinary clinics, laboratories as well as in cosmetic salons,

In the lighting fixture an optical system with a halogen bulb and absorption filter is used - this arrangement ensures uniform light intensity distribution over the entire illuminated area and absorbs a substantial part of thermal radiation, thus increasing the comfort of the doctor's work,

Specification:

Two-level light intensity regulation,

Bulb life 1000 h,

Light intensity at a distance of 1m - 30000 lx,

Color temperature 4000° K,

Range of light intensity adjustment - 20%,

Color rendering index Ra 90

Size of the light spot at a distance of 1m - 120 mm,

Height adjustment of the light source in the range 115 - 170 cm,

Dimensions:

Base diameter: 60 cm,

Height: 190 cm,

Weight: 9.2 KG,

Has a current inspection and is ready for operation,

Technical Passport issued (Service Report) valid for 12 months,

Warranty:

6 months for the domestic market (Poland),

3 months for the international market,

Possibility to extend warranty for a surcharge to 12 months for the domestic and international market,

Warranty period does not cover bulbs and other light sources installed in the device,

Financing available (Poland only): Installments, Leasing, Loan,

