



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



PATIENT HEATER MISTRAL-AIR PLUS MA1100-EU (Reconditioned)

Reconditioned (used),

Technical condition: very good,

Visual condition: very good,

Real product photos,

Dutch manufacture,

Power supply: 230 V,

Frequency: 50/60 Hz,

Power: 950 W,

PATIENT HEATER MISTRAL-AIR PLUS

combines modern design with power and simplicity. Provides quiet, effective, safe and reliable heating for your patients. Default initial temperature 38°C guarantees immediate start of therapy and reduces the risk of human error. In addition to overtemperature alarm, the too-low-temperature alarm is a standard built-in safety coefficient guaranteeing delivery of the set temperature to the patient. The standard patient hose fits blankets from all major manufacturers, both disposable and reusable,

Specification:

Forced air circulation heating,

Starting temperature: 38°C,

Set temperature: 32°C, 38°C or 43°C,

High and low temperature alarm function,

Filtration (HEPA): 0.3µm 99.99% / 0.2µm 99.86%,

Device mounted on a mobile stand with basket,

Five-spoke wheeled base,

Included:

English language user manual (PDF),

Power cord,

Dimensions:

Overall: 55 x 55 x 100 cm,

Hose length: 1.5 m,

Weight: 12.5 KG

Has current inspection and is ready for use,

Issued Technical Passport (Service Report) valid for 12 months,

Warranty:

6 months for domestic market (Poland),

3 months for international market,

Possibility to extend warranty for an additional fee up to 12 months for domestic and international markets,

Financing options (Poland only): Installments, Leasing, Loan,

